


SMALL CHILD CUFF (XS) FOR MOBI

Ref. 513710



14/06/2025 00:24:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

SMALL CHILD CUFF (XS) FOR MOBI



Ref. 513710

Description

Sizes available:

- XS arms with a circumference of 9 to 16 cm (513710)
- S arms from 15 to 26 cm in circumference (513711)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446027995	HS Code	90189010
Warranty (months)	12	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.143
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1	US Size (XXS - XXL)	XS

14/06/2025 00:24:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr