

PURPLE SPENGLER (ML) CUFF (22-42CM) FOR AUTOTENSIO

Ref. 527343



13/06/2025 23:46:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

PURPLE SPENGLER (ML) CUFF (22-42CM) FOR AUTOTENSIO



Ref. 527343

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446007522	HS Code	90262040
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Packed unit / Height (cm)	3
Packed unit / Length (cm)	16	Packed unit / Weight	0.124
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


 13/06/2025 23:46:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr