

CONICAL LOWER LIMB CUFF WITH 1 BLADDER FOR DDM TOURNIQUET

Ref. GCCS106

DDM

CE



 13/06/2025 22:52:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

CONICAL LOWER LIMB CUFF WITH 1 BLADDER FOR DDM TOURNIQUET



Ref. GCCS106

Those cuffs are autoclavable!

Description

Cuffs and lower limb cuff for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of **1 removable silicon bladder** for easy autoclaving.

Those ones are designed in a **conical form**

Several sizes are available

Dimensions of the bladder in cm

Lower limb cuff adult (GCCS104)	69 x 9 cm
Lower limb cuff adult L (GCCS105)	81 x 9 cm
Lower limb cuff adult XL (GCCS106)	85,5 x 12 cm

Dimensions de la poche en cm

Cuissard adulte (GCCS104)	69 x 9 cm
Cuissard adulte L (GCCS105)	81 x 9 cm
Cuissard adulte XL (GCCS106)	85,5 x 12 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203575	HS Code	90189084
Warranty (months)	24	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

13/06/2025 22:52:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr