

=====

Ref. 337849



 14/06/2025 00:17:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr



=====

Ref. 337849

Description

Single tube Clinic® cuff for Lian® + clamp compatible w/ Spengler hand-held BPMs Protocol patient monitors

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446030865	HS Code	9018901000
Warranty (months)	6	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	17	Packed unit / Height (cm)	8
Packed unit / Length (cm)	24	Packed unit / Weight	0.25
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

14/06/2025 00:17:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr