

MONOBLOC CONICAL ADULT ARM CUFF WITH ONE BLADDER AND A BUCKLE FOR DDM TOURNIQUETS

Ref. GBCMS105

DDM


CE



 13/06/2025 23:01:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

MONOBLOC CONICAL ADULT ARM CUFF WITH ONE BLADDER AND A BUCKLE FOR DDM TOURNIQUETS



Ref. **GBCMS105**

The monobloc conical arm cuff with one bladder and a buckle !

Description

This **Adult L** monobloc arm cuff made of **PU** is decontaminable. With a **conical design** and only **one bladder** it is made of a monobloc unit with a **buckle** to secure the attachment.

Width of the bladder : 10cm

Length of the bladder : 47cm

Total length : 70cm

The arm cuff can not be sterilized. It can be dry-cleaned or immersed with paying attention to maintain the connector of the tube out of water.

It is also available as lower limb cuffs

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468208129	HS Code	90189084
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			

13/06/2025 23:01:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr