

## MONOBLOCK STRAIGHT ARM CUFF WITH 1 BLADDER AND FASTENING STRAP FOR DDM PNEUMATIC TOURNIQUETS

Ref. GBDMS101

DDM

CE



14/06/2025 00:05:44 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## MONOBLOCK STRAIGHT ARM CUFF WITH 1 BLADDER AND FASTENING STRAP FOR DDM PNEUMATIC TOURNIQUETS



Ref. **GBDMS101**

Monoblock straight arm cuff with one bladder and fastening strap!

### Description


Monobloc PU decontaminable arm cuff are available in straight design. Made of a single monobloc unit with one bladder and a fastening strap to secure the attachment.



They are available in a range of different sizes easily indentifiable by color coding of the tubing. (A : bladder width, B: bladder length, C: total length)

		A	B	C
• New born arm cuff	White tube	3.5	29	38
• Child arm cuff	Pink tube	4.5	35	48
• Small adult arm cuff	Parme tube	6	46	61
• Adult arm cuff	Sky-blue tube	8	46	67
• Large adult arm cuff	Turquoise tube	10	47	70

Monobloc arm cuffs can not be sterilized. They can be dry cleaned or immersed paying attention to maintain the tube connector out of water.

 14/06/2025 00:05:44 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## MONOBLOCK STRAIGHT ARM CUFF WITH 1 BLADDER AND FASTENING STRAP FOR DDM PNEUMATIC TOURNIQUETS





Ref. **GBDMS101**

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700468208877	HS Code	90189084
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40
OTHER FEATURES			

 14/06/2025 00:05:44 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)