

## MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

medi Plogs


Ref. WSAM\_35BF



 13/06/2025 23:11:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

# MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE



Ref. **WSAM\_35BF**




## Description



- Clogs suitable for modern hospitals, operating theatres and intensive care units.
- Products designed and manufactured in collaboration with researchers and doctors.
- Antistatic clogs with non-slip and non-chafing soles
  - The back of the heel is sloped and rounded to make walking easier
  - Ventilation system allows air to circulate
  - Toe prints to enhance hold and spurs to aid the circulation
  - 100% CFC-free polyurethane
  - Washable by hand or machine at 60°C
  - Sizes: 35 - 47

## Features

CE Marking	NONE	Medical device	No
Bar code	3661040099705	HS Code	64029993
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	22	Packed unit / Height (cm)	12
Packed unit / Length (cm)	31	Packed unit / Weight	0.7
Bare Unit / Weight (kg)	1		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Navy blue	Shoe Size	35

 13/06/2025 23:11:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)