

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

medi Plogs

Ref. WSAM_35L



 13/06/2025 23:31:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE



Ref. **WSAM_35L**

=====

Description

Clogs suitable for modern hospitals, operating theatres and intensive care units.

Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

Features

CE Marking	NONE	Medical device	No
Bar code	3661040099712	HS Code	64029993
Country	Germany		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.16
Bare Unit / Weight (kg)	1		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	fuchsia	Shoe Size	35
--------	---------	-----------	----

📅 13/06/2025 23:31:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31
✉ info@gsh-med.fr