

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

medi Plogs

Ref. WSAM_39V



📅 13/06/2025 23:19:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE



Ref. **WSAM_39V**

=====

Description

Clogs suitable for modern hospitals, operating theatres and intensive care units.

Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

Features

| | | | |
|------------|---------------|----------------|----------|
| CE Marking | NONE | Medical device | No |
| Bar code | 3661040099972 | HS Code | 64029993 |
| Country | Germany | | |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|----|---------------------------|------|
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 1.16 |
| Bare Unit / Weight (kg) | 1 | | |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

| | | | |
|--------|-------|-----------|----|
| Colors | Green | Shoe Size | 39 |
|--------|-------|-----------|----|

📅 13/06/2025 23:19:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31
✉ info@gsh-med.fr