

LUXASCOPE OPHTHALMOSCOPE LED 2.5 V - LUXAMED

LUXAMED®

Ref. F1.416.114



25/06/2025 00:14:23 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LUXASCOPE OPHTHALMOSCOPE LED 2.5 V - LUXAMED

LUXAMED®

Ref. F1.416.114

Direct ophthalmoscope in pocket size

Description

- Direct ophthalmoscope in pocket size
- Refraction balance: + 20 bis - 20 D (± 1 2 3 4 6 8 10 15 20 D)
- 5 different apertures: reticle, large circle, small circle, semicircle, green filter
- Illumination: LED
- Can be dimmed for an individual glare-free brightness
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 50.000 hours by constant current drive
- Colour temperature: approx. 4.200 K
- Optimized design of head
- Two-component handle (aluminium/fibre-glass reinforced), also combinable with LuxaScope LED otoscope- and dermatoscope head
- Electronic ON/OFF push-button switch
- Power supply: 2x AAA Alkaline batteries, inclusive

- Available in :
- White : C1.516.914
 - Blue : C1.516.934
 - Pink : C1.516.974
 - Yellow : C1.516.944
 - Kids : C1.916.914
 - Black : C1.516.114
 - Green : C1.516.954
 - Purple : C1.516.964

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381964315	HS Code	90185090
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr