

## LUXASCOPE OPHTHALMOSCOPE LED 2.5 V - LUXAMED

**LUXAMED®**

Ref. F1.416.934



 24/06/2025 23:51:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

# LUXASCOPE OPHTHALMOSCOPE LED 2.5 V - LUXAMED

LUXAMED®

Ref. F1.416.934

Direct ophthalmoscope in pocket size

## Description

- Direct ophthalmoscope in pocket size
- Refraction balance: + 20 bis - 20 D (± 1 2 3 4 6 8 10 15 20 D)
- 5 different apertures: reticle, large circle, small circle, semicircle, green filter
- Illumination: LED
- Can be dimmed for an individual glare-free brightness
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 50.000 hours by constant current drive
- Colour temperature: approx. 4.200 K
- Optimized design of head
- Two-component handle (aluminium/fibre-glass reinforced), also combinable with LuxaScope LED otoscope- and dermatoscope head
- Electronic ON/OFF push-button switch
- Power supply: 2x AAA Alkaline batteries, inclusive

- Available in :
- White : C1.516.914
  - Blue : C1.516.934
  - Pink : C1.516.974
  - Yellow : C1.516.944
  - Kids : C1.916.914
  - Black : C1.516.114
  - Green : C1.516.954
  - Purple : C1.516.964

## Features

|                             |                |                           |          |
|-----------------------------|----------------|---------------------------|----------|
| CE Marking                  | CE             | Medical device            | Class I  |
| Bar code                    | 04260381964384 | HS Code                   | 90185090 |
| Warranty (months)           | 24             | Country                   | Germany  |
| DIMENSIONS AND WEIGHT       |                |                           |          |
| Packed unit / Width (cm)    | nc             | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)   | nc             | Packed unit / Weight      | nc       |
| SPECIFIC STORAGE CONDITIONS |                |                           |          |

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)