


## CONVERTIBLE BACK PACK EIFFEL, DE BOISSY

Ref. MDEIFF03

*De Boissy*  
— DEPUIS 1982



24/06/2025 23:56:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CONVERTIBLE BACK PACK EIFFEL, DE BOISSY

*De Boissy*  
— DEPUIS 1982

Ref. MDEIFF03

**It is worn in a backpack, hand or shoulder; A real Swiss army knife!**

### Description

Of a mottled gray, the polyester fabric is resistant, durable, water-repellent and therefore washable, it is accompanied by a reflective strip for safety, and YKK® slides for their quality. - Its main space

is modular and contains 4 zipped compartments, extractable (fixed with velcro): 3 transparencies to quickly identify the contents, and 1 isothermal, supplied with its reusable cold gel pocket, to slip up to 20 bulbs

-An extension zip makes it possible to enlarge this main compartment to store more equipment

- A second padded space can store computer / tablet and documents

- This model also has the large ventral pocket (45x21 cm) which opens completely with on one side 2 zipped storage nets and 2 flat pockets, and on the other side several elastic loops that allow to store small care supplies: pliers, scissors, pen lamp, thermometer, bandages, pens, etc ...

-Two external zipped pockets, plus one with elastic dedicated to the bottle of hydroalcoholic gel

-The shoulder straps, handles and shoulder strap are padded, solid and concealable to adopt the wear that suits best

-Dimensions w x h x d: 30 x 45 x 15/19 cm

-Weight: 2.1 kg

-Maximum recommended content: 6 kg / 19-25 liters

- Sold empty

### Features

CE Marking	NONE	Medical device	No
Bar code	3661040218953	HS Code	4202129000
Warranty (months)	12	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	50	Bare Unit / Width (cm)	15
Packed unit / Height (cm)	35	Bare Unit / Height (cm)	45
Packed unit / Length (cm)	50	Bare Unit / Length (cm)	30
Packed unit / Weight	2.1	Bare Unit / Weight (kg)	2.1
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

24/06/2025 23:56:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)