

## MOBI WITH ADULT ARMBAND (M) SPENGLER, 4 COLORS TO CHOOSE FROM



Ref. 513703



CE  
0459



 08/07/2025 23:11:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## MOBI WITH ADULT ARMBAND (M) SPENGLER, 4 COLORS TO CHOOSE FROM



Ref. 513703

Elegant and colorful, the MOBI blood pressure monitor revolutionizes the aneroid blood pressure monitor thanks to its new patented connectors.

### Description

The "LIFE LINK®" tubing allows the cuff to be changed in less than a second thanks to its magnetic connectors that ensure the connection and sealing between the pressure gauge and the tubing.

- AMBIDEXTRE : Mobi is specially designed to be ambidextrous, choose your configuration and change the position of the decompression valve, the "LIFE LINK®" connectors and the spoon.

- ERGONOMIC PEAR: The ergonomic pear is suitable for all hand morphologies and offers a volume gain of 11% on traditional pears.

- PRECISION VALVE: MOBI is equipped with a decompression valve whose shape is studied to offer the best sensitivity during decompression but also a very good ergonomics and grip.

- SHOCK RESISTANT: Due to its ABS plastic construction, its double inner shell and its elastomer rear part, MOBI is specially designed to offer great resistance to shocks and falls in everyday life.

- EASY ARMBAND with the revolutionary "LIFE LINK®" connectors, designed in TPU coating to be lighter and easy to clean.

- COLORED: Mobi exists in 4 colors: Carbon, Blueberry, Coral and Powder Pink-

Comes with an M cuff and a LIFE LINK® tubing

- Optional: 5 sizes of armbands (XS, S, M, L, XL)

### Features

|                             |               |                           |                               |
|-----------------------------|---------------|---------------------------|-------------------------------|
| CE Marking                  | CE0459        | Medical device            | Class Im (measuring function) |
| Bar code                    | 3700446027988 | HS Code                   | 90189100                      |
| Warranty (months)           | 36            | Country                   | China                         |
| DIMENSIONS AND WEIGHT       |               |                           |                               |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc                            |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 0.52                          |
| SPECIFIC STORAGE CONDITIONS |               |                           |                               |
| Minimum temperature         | -20           | Max temperature           | 70                            |
| Minimum humidity            | 20            | Maximum humidity          | 85                            |
| OTHER FEATURES              |               |                           |                               |
| Colors                      | Powder pink   |                           |                               |

08/07/2025 23:11:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)