

Ref. 513252



03/07/2025 22:13:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)



Ref. 513252

## Description

Efficiency and high-performance improved with the Clinic® version. This version of the Lian® Nano is lightweight and easy to handle and efficient thanks to its Clinic® cuff, suitable for liquid decontamination. Ideal for septic isolation or for preventing cross-contamination. Supplied in a nylon pouch.

## Features

|                             |               |                           |          |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking                  | CE            | Medical device            | No       |
| Bar code                    | 3700446024727 | HS Code                   | 90189010 |
| Warranty (months)           | 24            | Country                   | France   |
| DIMENSIONS AND WEIGHT       |               |                           |          |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 0.503    |
| SPECIFIC STORAGE CONDITIONS |               |                           |          |
| Minimum temperature         | -20           | Max temperature           | 70       |
| Minimum humidity            | 85            | Maximum humidity          | 85       |
| OTHER FEATURES              |               |                           |          |
| Packed unit / Quantity      | 1             |                           |          |

03/07/2025 22:13:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)