

F.O. MACINTOSH BLADE WITH INTEGRATED LIGHT GUIDE, LUXAMED ®,

Ref. E1.412.012

LUXAMED®



03/07/2025 22:26:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

F.O. MACINTOSH BLADE WITH INTEGRATED LIGHT GUIDE, LUXAMED ®,

LUXAMED®

Ref. E1.412.012

LUXAMED MACINTOSH laryngoscope blades with integrated fibre optics.

Description

Laryngoscope blade with integrated optical glass

- High-quality stainless steel construction for a long service life
- No dust collecting cracks – easy cleaning and sterilisation
- Maintenance-free
- 8.000 fibres to achieve the maximum service life and illumination
- Optimal light emission due to the largest possible cross section (3 x 7 mm)
- Brightness approx. 14.000 Lux at 3.7 V LED and 2.5 V LED
- Compatible with all handles according to ISO 7376 (green standard)
- Autoclavable up to 134°C – approx. 4.000 times
- Data-Matrix with most important information

Differents size available :

- 0 : E1.410.012
- 1 : E1.411.012
- 2 : E1.412.012
- 3 : E1.413.012
- 4 : E1.414.012
- 5 : E1.415.012

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960515	HS Code	90189084
Warranty (months)	24		

 **MENSIONS AND WEIGHT** *For information purposes only, can be subject to change*

Depth unit (Width) (mm)	Depth unit (Height) (mm)
GSH 30 rue Jean de Guiramand 13290 Aix-en-Provence France  www.gsh-med.fr  +33(0)4 42 90 31 31  info@gsh-med.fr	
Size 0	Size 2