

REFLEX HAMMER ACCORDING TO TROMNER, ALUMINIUM - LUXAMED

LUXAMED®

Ref. M1.111.915



03/07/2025 23:12:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

REFLEX HAMMER ACCORDING TO TROMNER, ALUMINIUM - LUXAMED

LUXAMED®

Ref. M1.111.915

Reflex Hammer according to Tromner | large

Description

- Reflex hammer made of aluminium / brass
 - Optimal weight distribution
 - With neurological sensor tip
 - Length: 23 cm | weight: 245 g
 - Without brush due to hygienic reasons
- Available in white (M1.111.915) and black (M1.111.112)

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960317	HS Code	90189084
Warranty (months)	24		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	White		

📅 03/07/2025 23:12:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr