

REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

LUXAMED®

Ref. M1.131.914



 03/07/2025 23:25:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

LUXAMED®

Ref. M1.131.914

Reflex Hammer according to Buck | small

Description

- Reflex hammer made of aluminium
- Optimal weight distribution
- Powder coating
- With neurological sensor tip
- Length: 18 cm | weight: 90 g
- Without brush due to hygienic reasons

Available in :

White: M1.121.914

Blue :M1.121.812

Pink :M1.121.612

Yellow: M1.121.513

Kid : M1.131.914

Black M1.121.112

Green: M1.121.412

Purple: M1.121.312

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381963356	HS Code	90189084
Warranty (months)	24		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	40
Minimum humidity	40	Max humidity	60

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr