


## CLAR HEAD MIRROR-ADJUSTABLE HEADBAND FOR ALL SIZES+POUCH



Ref. 570150



10/07/2025 00:12:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CLAR HEAD MIRROR-ADJUSTABLE HEADBAND FOR ALL SIZES+POUCH



Ref. 570150

### Description

Suitable for all head sizes thanks to a flexible metal headband.  
 Indirect lighting by reflecting mirror.  
 Suited to lighting required in daily practice.  
 Built-in battery compartment, two CR 123 batteries  
 Supplied with carrying bag

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446011390	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	23	Packed unit / Height (cm)	11
Packed unit / Length (cm)	24	Packed unit / Weight	0.655
Bare Unit / Weight (kg)	0.23		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

10/07/2025 00:12:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
 13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)