


IDEAL+ SINGLE HEAD ADULT,BLUE

Ref. DH7081_R



15/07/2025 21:47:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

IDEAL+ SINGLE HEAD ADULT,BLUE



Ref. DH7081_R

Description

- Stethoscope with single adult metal horn (O: 45 mm)•
- Anti-cold ring and horn of the color of the tube•
- Aluminum lyre incorporated into the tubing
- Also available in double pavilion for adults and a pediatric model

Available in Black (DH7081_N) or Red (DH7081_R)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041000137	HS Code	9018908400
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.15
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Red	Packed unit / Quantity	1

15/07/2025 21:47:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr