

BOX OF 800 PROBE COVERS FOR TEMPO EAR PRO AND TEMPO EAR PRO +



Ref. THE003ACC1



17/07/2025 21:03:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

BOX OF 800 PROBE COVERS FOR TEMPO EAR PRO AND TEMPO EAR PRO +



Ref. THE003ACC1

- Box of 800 probe covers for Tempo EAR PRO and Tempo EAR PRO +

Description

- Box of 800 probe covers for Tempo EAR PRO and Tempo EAR PRO +
- Dispenser function
- Pack of 20

Features

CE Marking	CE	Medical device	Class I
Bar code	6942284605264	HS Code	3926909090
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-25	Max temperature	55
Minimum humidity	15	Maximum humidity	95

OTHER FEATURES

17/07/2025 21:03:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr