

Ref. 334121



 24/07/2025 21:07:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr



Ref. 334121

Description

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Metal, Lian® Nano, Lian® Scenic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446001995	HS Code	90189010
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.11
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

24/07/2025 21:07:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr