

Ref. 334150



24/07/2025 21:19:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)



Ref. 334150

## Description

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Metal, Lian® Nano, Lian® Scenic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

## Features

|                             |               |                           |            |
|-----------------------------|---------------|---------------------------|------------|
| CE Marking                  | CE            | Medical device            | Class I    |
| Bar code                    | 3700446026622 | HS Code                   | 9018901000 |
| Warranty (months)           | 6             | Country                   | Indonesia  |
| DIMENSIONS AND WEIGHT       |               |                           |            |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc         |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 0.335      |
| SPECIFIC STORAGE CONDITIONS |               |                           |            |
| Minimum temperature         | NA            | Max temperature           | NA         |
| Minimum humidity            | NA            | Maximum humidity          | NA         |
| OTHER FEATURES              |               |                           |            |
| Packed unit / Quantity      | 1             |                           |            |

24/07/2025 21:19:39 - For information purposes only, can be subject to change