

IDEAL+ DOUBLE HEAD ADULT, BLUE

Ref. DH7073_RC



24/07/2025 21:16:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

IDEAL+ DOUBLE HEAD ADULT, BLUE



Ref. DH7073_RC

=====

Description

- Stethoscope with double metal horn adult (O: 45/30 mm)•
- Anti-cold ring and horn of the color of the tube•
- Aluminum lyre incorporated into the tubing
- Also available in single pavilion for adults

Available in Blue (DH7073_B), Black (DH7073_N), Red (DH7073_R) or Pink (DH7073_RC) for the paediatric version (O: 35/30 mm)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041000106	HS Code	9018908400
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Packed unit / Height (cm)	4
Packed unit / Length (cm)	25	Packed unit / Weight	0.15
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Pink	Packed unit / Quantity	1

24/07/2025 21:16:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr