

PERFORATED BOXES, 18 X 8 X 4 CM, PERFORATED TOP / BOTTOM

Ref. AI18084P



 30/07/2025 09:32:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

PERFORATED BOXES, 18 X 8 X 4 CM, PERFORATED TOP / BOTTOM



Ref. AI18084P

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067728	HS Code	90189084
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	21	Bare Unit / Width (cm)	8
Packed unit / Height (cm)	10	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	23	Bare Unit / Length (cm)	18
Packed unit / Weight	0.248	Bare Unit / Weight (kg)	0.248
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 09:32:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr