

## S/S BOX PERFORATED 20 X 10 X 3 CM PERFORATED TOP / BOTTOM

Ref. AI20103P



30/07/2025 03:17:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## S/S BOX PERFORATED 20 X 10 X 3 CM PERFORATED TOP / BOTTOM



Ref. AI20103P

### Features

|                             |               |                         |          |
|-----------------------------|---------------|-------------------------|----------|
| CE Marking                  | CE            | Medical device          | Class I  |
| Bar code                    | 3661040067735 | HS Code                 | 90189084 |
| Warranty (months)           | 12            | Country                 | Pakistan |
| DIMENSIONS AND WEIGHT       |               |                         |          |
| Packed unit / Width (cm)    | 13            | Bare Unit / Width (cm)  | 13.6     |
| Packed unit / Height (cm)   | 5             | Bare Unit / Height (cm) | 5.5      |
| Packed unit / Length (cm)   | 28            | Bare Unit / Length (cm) | 28       |
| Packed unit / Weight        | 0.325         | Bare Unit / Weight (kg) | 0.325    |
| SPECIFIC STORAGE CONDITIONS |               |                         |          |
| Minimum temperature         | NA            | Max temperature         | NA       |
| Minimum humidity            | NA            | Maximum humidity        | NA       |
| OTHER FEATURES              |               |                         |          |
| Packed unit / Quantity      | 1             |                         |          |

30/07/2025 03:17:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)