

STAINLESS STEEL BOX, 25 X 10 X 5 CM - PERFORATED TOP/BOTTOM

Ref. AI25105P



30/07/2025 01:38:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

STAINLESS STEEL BOX, 25 X 10 X 5 CM - PERFORATED TOP/BOTTOM



Ref. AI25105P

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067766	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	26	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	12	Bare Unit / Length (cm)	25
Packed unit / Weight	6.7		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 01:38:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr