

S/S BOX PERFORATED 25 X 12 X 6 CM PERFORATED TOP / BOTTOM

Ref. AI25126P



29/07/2025 23:33:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

S/S BOX PERFORATED 25 X 12 X 6 CM PERFORATED TOP / BOTTOM



Ref. AI25126P

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067773	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Bare Unit / Width (cm)	12.9
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	6.6
Packed unit / Length (cm)	26	Bare Unit / Length (cm)	25.7
Packed unit / Weight	0.67	Bare Unit / Weight (kg)	0.613
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 23:33:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr