

S/S BOX 30 X 12 X 6 CM

Ref. AI30126

 **Holtex**



CE



 30/07/2025 01:42:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

S/S BOX 30 X 12 X 6 CM



Ref. AI30126

Convenient for storing surgical instruments

Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019598	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Bare Unit / Width (cm)	12
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	6
Packed unit / Length (cm)	30	Bare Unit / Length (cm)	30
Packed unit / Weight	0.81	Bare Unit / Weight (kg)	0.758
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 01:42:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr