

S/S BOX PERFORATED 35 X 18 X 8 CM PERFORATED TOP / BOTTOM

Ref. AI35188P



29/07/2025 23:49:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

S/S BOX PERFORATED 35 X 18 X 8 CM PERFORATED TOP / BOTTOM



Ref. AI35188P

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067872	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	19	Bare Unit / Width (cm)	18
Packed unit / Height (cm)	9	Bare Unit / Height (cm)	8
Packed unit / Length (cm)	38	Bare Unit / Length (cm)	36
Packed unit / Weight	1.142	Bare Unit / Weight (kg)	1.142
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 23:49:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr