

S/S BOX PERFORATED 40 X 20 X 9 CM PERFORATED TOP / BOTTOM

Ref. AI40209P



 30/07/2025 00:24:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

S/S BOX PERFORATED 40 X 20 X 9 CM PERFORATED TOP / BOTTOM



Ref. AI40209P

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067889	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	21	Packed unit / Height (cm)	11
Packed unit / Length (cm)	43	Packed unit / Weight	1.51
Bare Unit / Weight (kg)	1.4		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 00:24:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr