

B/B SCISSORS STRAIGHT 24 CM

Ref. IC11424



30/07/2025 00:42:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

B/B SCISSORS STRAIGHT 24 CM



Ref. IC11424

Ideal for small care (dressings, compresses, etc.)

Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations.

Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016).

- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" leaflet

Features

| | | | |
|-----------------------------|---------------|-------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040083346 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Bare Unit / Width (cm) | 5.5 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 0.5 |
| Packed unit / Length (cm) | nc | Bare Unit / Length (cm) | 24.6 |
| Packed unit / Weight | 0.126 | Bare Unit / Weight (kg) | 0.118 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

30/07/2025 00:42:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr