

ASKIR C30 ASKIR C30 SUCTION UNIT UU FLOVAC 2X2L



Ref. BCASPI37



30/07/2025 10:34:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

ASKIR C30 ASKIR C30 SUCTION UNIT UU FLOVAC 2X2L



Ref. BCASPI37

Surgical Suction Unit on wheels and mobile, intended to drain body fluids and to practice oral, nasal and tracheal aspiration in adults and children.

Description

Ca-mi Askir C30 Suction Unit, Flovac® single-use pouches, 2 x 2000 ml

Surgical Suction Unit on wheels and mobile, intended to drain body fluids and to practice oral, nasal and tracheal aspiration in adults and children.

Intended for professional use and equipped with a powerful suction system that can drain up to 40 l/min max, Askir C30 is equipped with two collection jars and a five-caster system making it easy to transport.

Available in Basic and FS

versions- BASIC FLOVAC version, including: Fluid collection jar with anti-reflux valve / 8 x 14 Silicone hoses (autoclavable) and conical connectors o 10-11-12 mm / Antibacterial & hydrophobic filter (single patient) / Power cord with Schuko plug

-

Flovac® single-use pockets-

Main

system equipped with a pressure gauge (bar & kPa) and a vacuum regulator

- Engine: Oil-free piston pump requiring no maintenance

-

Power supply: 230V-50Hz

-

Maximum vacuum (adjustable): -0.80 Bar -80 kPa -600 mmHg

-

Max free air flow: 40 l/min

-

Noise level: 61 dBA

- Weight: 6.2 kg / Dimensions: 32 x 99 x 30 cm

Features



30/07/2025 10:34:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

Packed unit / Length (cm)	nc	Packed unit / Weight	9.4
Bare Unit / Weight (kg)	6.2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			