

X-RAY VIEWER STANDARD-SCREEN 45W

Ref. OTA1001



30/07/2025 02:22:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guirmand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

X-RAY VIEWER STANDARD-SCREEN 45W



Ref. OTA1001

Description

- High resistance aluminium alloy structure. Corners are protected by a steel plate covered in ABS Plexiglas polymer screen. High resistance aluminium alloy body. The interior and exterior are protected by a sprayed anti-static powder.

- These x-ray viewers have a radiograph fixing system that allows them to be placed and removed easily and without risk of damage. They can be free-standing or fixed to a wall

- Lamp Spec. T8

- 1 screen, size: 36 x 43.5 cm
- External dimensions: 44 x 58 x 12 cm

- Package size: 70 x 53 x 20 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040195858	HS Code	90105010
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	50	Bare Unit / Width (cm)	44
Packed unit / Height (cm)	15	Bare Unit / Height (cm)	10
Packed unit / Length (cm)	61	Bare Unit / Length (cm)	56
Packed unit / Weight	7.4	Bare Unit / Weight (kg)	6
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 02:22:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr