


EXTRA-FLAT X-RAY VIEWER - 2 PANELS

Ref. OTA1102



30/07/2025 04:04:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

EXTRA-FLAT X-RAY VIEWER - 2 PANELS



Ref. OTA1102

Description

These negatoscopes offer a mounting system for radios allowing simple installation and removal, and without consequence for radios. - The structure of the screen is made of high-strength aluminum alloy

- 2.4 cm thick - The angles are protected by a steel plate, covered with ABS
- Methacrylic polymer screen- The body of the negatoscope is made of high strength aluminum alloy
- The inside and outside are covered with an antistatic sprayed powder
- Voltage / Frequency: 220V / 30-100Hz
- Average brightness> 1500 cd / m2; >5,000 lux.
- Lighting life: 100,000 h
- Color temperature: 6500k- Screen dimensions: 74 x 43.5 cm
- External dimensions: 85 x 54.5 x 2.4 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040195896	HS Code	90105010
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	65	Bare Unit / Width (cm)	85
Packed unit / Height (cm)	10	Bare Unit / Height (cm)	2.4
Packed unit / Length (cm)	93	Bare Unit / Length (cm)	54.5
Packed unit / Weight	15.1	Bare Unit / Weight (kg)	13.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 04:04:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr