

S/S THERMOMETER JAR H.100 X 30 (DIAM) MM

Ref. AI11611



30/07/2025 06:01:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

S/S THERMOMETER JAR H.100 X 30 (DIAM) MM



Ref. AI11611

Description

Autoclavable stainless steel thermometer holder for medical use.

Medical device that allows a conventional medical thermometer to be temporarily stored without the risk of it falling and damage, and also avoids direct contact between the thermometer and a work surface to reduce the risk of cross-contamination.

Clean and disinfect well after each use.

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations.

Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019291	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	8	Bare Unit / Width (cm)	3
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	10
Packed unit / Length (cm)	11	Bare Unit / Length (cm)	10
Packed unit / Weight	0.088	Bare Unit / Weight (kg)	0.088
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 06:01:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr