

S. S BED PAN TYPE 37 X 29CM

Ref. AI13729



29/07/2025 22:55:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

S. S BED PAN TYPE 37 X 29CM



Ref. AI13729

Description

Capacity 5 Liters

Maximum weight supported: 90 kg

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019345	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	38	Bare Unit / Width (cm)	29
Packed unit / Height (cm)	32	Packed unit / Length (cm)	54
Bare Unit / Length (cm)	37	Packed unit / Weight	0.69
Bare Unit / Weight (kg)	0.68		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 22:55:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr