

S/S BED PAN WITHOUT COVER SLIPPER TYPE 48 X 30CM

Ref. AI14830



30/07/2025 03:07:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

S/S BED PAN WITHOUT COVER SLIPPER TYPE 48 X 30CM



Ref. AI14830

Description

Capacity 6 Liters

Maximum weight supported: 90kg

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019352	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	55	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	8	Packed unit / Length (cm)	30
Bare Unit / Length (cm)	48	Packed unit / Weight	1.08
Bare Unit / Weight (kg)	1.07		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 03:07:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr