

REPLACEMENT INFRARED EAR THERMOMETER FOR NC3/NC5

COMeN

Ref. BCTHER01



 30/07/2025 00:02:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

REPLACEMENT INFRARED EAR THERMOMETER FOR NC3/NC5

COMEN

Ref. BCTHER01

The wireless infrared thermometer that gives a measurement in 2 to 3 seconds!

Features

CE Marking	NONE	Medical device	No
Bar code	3661040219349	HS Code	9025199010
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.1		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 00:02:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr