

ANGULATED DISSECTING FORCEPS 16 CM

Ref. IP04516



30/07/2025 03:11:44 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

ANGULATED DISSECTING FORCEPS 16 CM



Ref. IP04516

Allows to carry out care and facilitate the gripping of small objects

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040021843	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.1
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.026	Bare Unit / Weight (kg)	0.025
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/07/2025 03:11:44 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr