

ARNEODO HOOK FOR PHACO EMULSIFICATION



Ref. IC60480



30/07/2025 04:05:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ARNEODO HOOK FOR PHACO EMULSIFICATION



Ref. IC60480

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040084046	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	1	Bare Unit / Width (cm)	1.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	12	Bare Unit / Length (cm)	12.1
Packed unit / Weight	0.016	Bare Unit / Weight (kg)	0.013
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 04:05:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr