

ANGLED ANTEROIR CHAMBER IRRIGATION CANNULA (0.9 X 0.35 MM)



Ref. IC92116



 30/07/2025 00:41:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

ANGLED ANTEROIR CHAMBER IRRIGATION CANNULA (0.9 X 0.35 MM)



Ref. IC92116

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040084114	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	4.4
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	5.8
Packed unit / Weight	0.014	Bare Unit / Weight (kg)	0.014
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/07/2025 00:41:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr