

DOYEN MOUTH-GAG ADULT, 13.5 CM

Ref. ID03602



29/07/2025 22:37:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DOYEN MOUTH-GAG ADULT, 13.5 CM



Ref. ID03602

Description

Dean stainless steel mouth opener for adults. Total length 13.5 cm (Hooks: 60 mm - Ends: 24 mm).

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040013459	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	14.3
Packed unit / Weight	0.121	Bare Unit / Weight (kg)	0.121
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 22:37:53 - For information purposes only, can be subject to change