

S/S KIDNEY DISH 30 X 14 X 4,5 CM

Ref. AI28124



29/07/2025 22:30:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

S/S KIDNEY DISH 30 X 14 X 4,5 CM



Ref. AI28124

Ideal for storing instruments or temporarily installing medical waste

Description

- Decontamination with detergent and disinfectant -

Autoclavable (134°-18 minutes)- In the "

ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual

- Reference 100% IN-STORE: this service instrumentation reference is guaranteed available immediately

- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019567	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	17	Bare Unit / Width (cm)	14
Packed unit / Height (cm)	7	Bare Unit / Height (cm)	4.5
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	30
Packed unit / Weight	0.28	Bare Unit / Weight (kg)	0.246
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 22:30:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr