

LEXER GOUGE 25CM 12MM

Ref. IG01135



29/07/2025 23:37:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LEXER GOUGE 25CM 12MM



Ref. IG01135

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040044040	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	29.5
Packed unit / Weight	0.151	Bare Unit / Weight (kg)	0.151
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 23:37:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr