

## MURPHY-PEAN FORCEPS STRAIGHT 14 CM

Ref. IP21414



30/07/2025 02:31:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## MURPHY-PEAN FORCEPS STRAIGHT 14 CM



Ref. IP21414

Used during surgical procedures (stop bleeding, hold compress, etc.)

### Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040023625	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	23	Bare Unit / Width (cm)	8.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.4
Packed unit / Length (cm)	9	Bare Unit / Length (cm)	14
Packed unit / Weight	0.04	Bare Unit / Weight (kg)	0.04
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 02:31:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)