

LOUTE LOOP TWISTER

Ref. IT10100



 30/07/2025 00:57:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LOUTE LOOP TWISTER



Ref. IT10100

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040031651	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Bare Unit / Width (cm)	7.8
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	4.3
Packed unit / Length (cm)	26	Bare Unit / Length (cm)	26
Packed unit / Weight	0.378	Bare Unit / Weight (kg)	0.378
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 00:57:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr