

VACHER SNARE

Ref. IS18431



30/07/2025 06:27:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

VACHER SNARE



Ref. IS18431

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040087764	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	18	Bare Unit / Width (cm)	18
Packed unit / Height (cm)	1	Bare Unit / Height (cm)	1
Packed unit / Length (cm)	35	Bare Unit / Length (cm)	35
Packed unit / Weight	0.2	Bare Unit / Weight (kg)	0.2
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 06:27:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr