

OCHSNER TROCART 18 FRENCH GAUGE (6 MM)

Ref. IT11618



30/07/2025 06:52:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

OCHSNER TROCART 18 FRENCH GAUGE (6 MM)



Ref. IT11618

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040031743	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	3	Packed unit / Height (cm)	2
Packed unit / Length (cm)	17	Packed unit / Weight	0.024
Bare Unit / Weight (kg)	0.024		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 06:52:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr