

## DOYEN ABDOMINAL RETRACTOR 40 MM X 85 MM

Ref. IV03485



29/07/2025 22:38:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DOYEN ABDOMINAL RETRACTOR 40 MM X 85 MM



Ref. IV03485

To expose the vaginal cavity

### Description

Holtex is specialist in stainless steel surgical instruments since 1980. Holtex instruments are all made of top quality stainless steel in order to guarantee good resistance to corrosion and durability.

> This reference is made with Stainless Steel AISI 304.

The Holtex surgical instruments are classified as medical devices and comply ISO 13485 standards and the 93/42/EC directive.

EN ISO 13485: 2016, ISO 14971: 2019, EN 1041: 2008, and EN ISO 15223-1

Holtex Instruments are guaranteed for 2 years under normal conditions of use. The warranty does not cover wear parts, shocks and maintenance defects made by the user

Packaging and traceability: Holtex instruments are supplied in individual bags, sealed and labelled with : Logo, and address of the manufacturer, the product reference, the detailed designation in French and English, the lot number ("LOT"), the CE marking, and the barcode.

Cleaning: Please read carefully instructions of cleaning and sterilisation of stainless-steel equipment inserted in each individual bag. Instructions can be also downloaded from the product page on our website.

### Features

|                   |               |                |          |
|-------------------|---------------|----------------|----------|
| CE Marking        | CE            | Medical device | Class I  |
| Bar code          | 3661040032320 | HS Code        | 90189084 |
| Warranty (months) | 24            | Country        | Pakistan |

#### DIMENSIONS AND WEIGHT

|                           |       |                         |       |
|---------------------------|-------|-------------------------|-------|
| Packed unit / Width (cm)  | nc    | Bare Unit / Width (cm)  | 8.6   |
| Packed unit / Height (cm) | nc    | Bare Unit / Height (cm) | 5     |
| Packed unit / Length (cm) | nc    | Bare Unit / Length (cm) | 27.5  |
| Packed unit / Weight      | 0.234 | Bare Unit / Weight (kg) | 0.234 |

#### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

#### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|

29/07/2025 22:38:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)