


ALUMINIUM TRAY COLOR 20 X 10 CM, AZUR

Ref. AAX201_B



29/07/2025 23:59:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ALUMINIUM TRAY COLOR 20 X 10 CM, AZUR



Ref. AAX201_B

Description

- Hot sterilization or autoclave (134°-18 minutes)
- Do not wash in an ultrasonic cleaner or in a washing machine
- Warning: aluminum is incompatible with some decontamination products

Format 20 x 10 x 2 cm

Available colors:

- Azure- Blue

-

Golden

- Grey

- Red

- Green

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067131	HS Code	90184990
Warranty (months)	12	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	10	Bare Unit / Width (cm)	10.2
Packed unit / Height (cm)	2	Bare Unit / Height (cm)	2.2
Packed unit / Length (cm)	20	Bare Unit / Length (cm)	20
Packed unit / Weight	0.074	Bare Unit / Weight (kg)	0.073
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Blue	Packed unit / Quantity	1

29/07/2025 23:59:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr